

When stones cause trouble

Endoscopical Large Stone Management in Gastroenterology and Urology

Lithotron EL 27 Compact State-of-the-Art Technology made in Germany Walz EHL Electrohydraulic Lithotripsy

.

Walz Lithotron EL 27 Compact

Large Stone Management


Where other equipment is limited our Walz EHL helps


WALZ

Ready to use: the complete Walz EL 27 Compact device

High Tech in touch with the stone

Fragmenting hard stones while treating soft tissue with care is like peeling raw eggs! And that's exactly what you can do with our Lithotron EL 27 Compact. Walz devices represent state-of-the-art-technology of the highest quality made in Germany.






before

Practical application of the Walz EL 27 in the common bile duct





The Technology:

Micro-spark plasmas are produced between bipolar electrodes at the centre of the tip of the flexible probes by short high voltage pulses. The rapidly expanding and collapsing spark plasma bubbles generate micro shock waves with steep edges and minimized tensile phases in the surrounding liquid. Additionally, micro liquid jets with high velocity are produced which are directed at the stone.

The Effects:

These micro shock waves in combination with the impacts of the micro liquid jets destroy the stone with highest dual effectiveness – while the minimized tensile phases and the short range of the micro shock waves help to treat the surrounding tissue with care.

40 years of experience in EHL technologies pay off!

EHL probes

- Can be ordered sterile or non-sterile
- Highly flexible but keep their position
- Length: 450 -3000 mm
- Diameter: 2 Fr., 2,4 Fr., 3 Fr., 4,5 Fr., 7 Fr., 10 Fr.
- Precisely worked probe tip for reproduceable shock waves
- Centred micro spark plasmas keep distances to the tissue and facilitate targeting
- Endoscope protection through rounded edges at the tip of the probe, no breaking problem compared to laser fiber

The Walz EL 27 Compact is the first choice when it comes to lithotripsy devices. Its features and benefits at a glance:

- **Highest efficiency** compared to other technologies, short operation times
- Stones are normally **disintegrated** with only a few shock waves
- **Security:** minimized tensile phases of the micro shock waves help to treat soft tissue with care.
- Highly localized energy transfer compared to laser and extracorporal shock waves
- Economical: maintenance-free device
- **Endoscope protection** through rounded probe tips; no breaking problem compared to laser fibers
- **Controlled energy levels:** 3 intensities (up to 950 mJ) adjusted to usage in gastroenterology and urology
- Very fast pulses, sharply rising amplitudes for **most** effective stone fragmentation

Technical Specifications:

Dimensions

Width: 320 mm, Height: 120 mm, Depth: 240 mm Weight: 7,5 kg

Power supply Voltage Versions	_ /-	
	115V <u>+</u> 10%	1,2 A
	230V ± 10%	0,6 A

EHL Electric charge

Intensity A/Pulse Frequency (max.) 250 mJ/60 Hz Intensity B/Pulse Frequency (max.) 500 mJ/50 Hz Intensity C/Pulse Frequency (max.) 950 mJ/40 Hz

WALZ ELEKTRONIK GMBH

- Located in Rohrdorf near Stuttgart/ Germany.
- Development, production and service of lithotripsy devices for Urology, Gastroenterology and Industry.
- Pioneers in electrohydraulic lithotripsy (EHL) with nearly 40 years of experience.
- · More than 20 property right (patent) applications.
- · More than 2.000 devices sold.
- Certifications
- · ISO 13485
- Appendix II of the directive 93/42/EWG

History

- **1973** Diplom-Ingenieur Volker Walz develops the first EHL device and related probes at the University of Stuttgart.
- **1974** Walz sells the first commercial EHL device.
- **1994** The first combined EHL/EKL (electrokinetic lithotripsy) device enters the market.
- **2003** The LithoRapid is an EKL-only device with improved effectivity.
- **2012** Bernd Vollmer enters Walz Elektronik GmbH as successor of Volker Walz.

Walz Elektronik GmbH Walddorfer Str. 40 72229 Rohrdorf Germany T +49 (0) 7452 2020 F +49 (0) 7452 3826 info@walz-el.de www.walz-el.de

